
BENMOR MEDICAL LAUNCHES NEW AURUM RANGE

BENMOR MEDICAL IS PROUD TO ANNOUNCE THE INTRODUCTION OF 'THE AURUM RANGE'.
NEW PRODUCTS INCLUDE THE AURUM EXPANDABLE BED



AURUM EXPANDABLE BED

The Aurum bariatric bed has been designed specifically for size and strength. The width adjustment feature allows the bed to be used alongside standard beds during everyday use but can be expanded quickly and simply to accommodate bariatric patients. The width adjustment also enables the bed to be transported as a standard sized bed therefore will fit through doors.

The four section, fully profiling mattress platform is powered by electric actuators which enables the bed to cope with patients up to 65st/414kg.

The integral scale unit enables staff and carers to carry out accurate patient weights up to 65st/414kg.

The bed has also been fitted with high strength castors and an integral linen tray.

