

bradshaw bariatric bed

bariatric community bed



Key features / spec

320kg (50st) swl

Platform height from 259mm -
625mm (10" - 25")

1170mm (46") wide, caters
for wider patients

4 section bed easy to
assemble/dismantle

Electronically operated with
trendelenburg functions

Easy to use handset controller

The bradshaw bariatric community bed has been designed for use in the community, nursing and residential homes. Compared to a hospital bed the community bed has a more domestic appearance due to the wooden head and footboards.

The bed has a three section, fully powered mattress platform and is ideal for carers dealing with bariatric patients in a carehome or their own residence.

Product code: -
BMPB050