

S-capepod[®]

evacuation device



The S-capepod provides a quick and safe evacuation of bed-bound and immobile patients in an emergency – where speed and ease of deployment is critical.

Features	Benefits
300kg (47.3st) max load	Able to cater for a range of patient weights
Fast and safe evacuation	You are in control of the evacuation of patients in an emergency
2 - person operation advised	Requires very little training
Complies with current standards	Ensure your evacuation plans are in place
Virtually no maintenance	Strong hard wearing and easy to clean

The S-capepod is a simple-to-deploy evacuation device, stored below the bed mattress, which when required can be simply unfolded and fitted around the person laid on the bed mattress.

Using the folding flap sections at each side and at the foot-end of the bed to create a sturdy bag in which to envelop the patient and mattress.

Hand-hold straps at the head and foot-end allow the mattress and occupant to be easily maneuvered along the floor and down stairs.

The S-capepod is suitable for loads up to 300kg (47st).

Product codes

BMPBSCAPE3
BMPBSCAPE5

Bariatric S-Capepod 300kg SWL
Bariatric S-Capepod 500kg SWL