

bed scale

portable bed scale with up to 1000kg weight capacity

The bed scale has been designed as a convenient solution to patient weighing, where the patient can only be weighed whilst still in bed, and where moving the bed to weigh is not practical.



- Portable scale for bedridden patients
- Preset tare for quick and easy patient weighing
- Low profile design compatible with almost any bed
- High capacity at 1000kg
- Highly accurate at 500g

Class III Approval	
Capacity	1000kg
Graduations	500g
Adaptor	3 pin tip + ve 12V 1A
Power Supply	Mains / Rechargeable



Product code: BMW005

Quick, easy and accurate bedridden patient weighing for hospitals and care homes

Being ultra-low profile, the four pads of the bed scale are small in size, making this our most lightweight bed weighing solution yet.

Secondly, the connectors for each pad are colour coded, meaning setting up the scale for weighing a patient in their bed is quick and easy.

The scale provides an accurate reading of the patient's weight, without the weight of the bed registered. It achieves this via the Pre-Set Tare; once the indicator knows the weight of the beds you use, it will remember it each and every time you use the scale thereafter.

This scale allows a patient to be accurately weighed whilst still in their bed, using four pads that are placed under the bed's wheels. The pads are ultra low profile, meaning they are suitable for almost any bed.

Each of the four pads, weighing at just 5kg each, can be transported between rooms, wards or floors with the use of the carry cases, or the optional trolley. There are also carry handles on each of the four pads for extra portability.



Product Code
BMW005

Available for purchase only



ISO 9001:2015 - FM 607485
ISO 13485:2003 &
ISO EN 13485:2012 - MD 607484