

high capacity chair scale

ideal for bariatric patients, hospital wards & care homes

The chair scale with a high capacity and super-sized seat, suitable for weighing bariatric patients. The chair scale has hinged armrests and footrests for patient comfort and easy patient transfer. Two rear brake wheels and fixed wheels at the front allow for smooth movement when a patient is being transferred.



- Extra Wide Super Sized Seat
- Robust Construction
- Suitable for Patients of all Sizes and Weights
- Ideal for Bariatric Patients
- Optional Handrail
- Optional Footrests



Product code: BMW007

Capacity	300kg
Graduations	50g<150kg>100g
Seat Dimensions	530mm x 460mm
Weight of scale	27kg
Power Supply	Mains / Rechargeable