

dynamic air cushion

pressure relieving bariatric air cushion



Key features / spec

Two way stretch material, conforms to the contours of the patient's body

Various widths available, compatible with our seating range

Adjustable pressure, provides optimum pressure relief for patients

24 hours air retention when pump is disconnected

Waterproof, wipe clean cover for easy cleaning and infection control

The dynamic air cushion has been designed to give high risk pressure relieving care to bariatric patients who suffer from pressure sores and general discomfort.

The cushions have an alternation cycle of ABAB and operate using a pump with adjustable pressure. The cover is made from a two way stretch material which conforms to the contours of the body and is vapour permeable, water resistant, fire retardant and soft whilst still being durable.

Product codes

BMAC001-5M - 25" 349kg (55st) swl

BMAC002-5M - 27" 349kg (55st) swl

BMAC003-5M - 29" 349kg (55st) swl

BMAC004-5M - 31" 349kg (55st) swl