



# Evacuation Equipment Guide

## Solutions for safe evacuation

# About Us

**Benmor Medical** provide the highest quality equipment and services to the healthcare industry. Unique to the market, we design and manufacture innovative products which help to improve patient care and maintain safety for healthcare professionals.

We were the first company in the UK to specialise in bariatric patient handling equipment.

In the UK we provide a Nationwide 24/7 service providing emergency rentals to hospitals and the community.

In addition to delivery and installation, our highly skilled technicians provide full training on all of our equipment.

“*Benmor Medical provide a fast and efficient service that enabled us to care for patients in a bed most suitable for their needs*” - Moving & Handling Co-ordinator

## Benmor Medical Delivers:



24/7 Emergency Nationwide Rental Equipment



Expert Advice & Innovative Products



In House Clinical Support



Patient Assessments & Weighing Service



Training & Installation



Products Designed & Manufactured in the UK



CPD Certified Training Courses



Key Focus on Compliance & Ethics



0333 800 9000



[info@benmormedical.co.uk](mailto:info@benmormedical.co.uk)

# Evacuation Equipment Range



**TETCON**

In an emergency situation, it is crucial that everyone in the building can reach a point of safety as quickly and safely as possible.

With this in mind, Netherlands based textile company Tetcon have developed an extensive range of evacuation equipment.

Equipment which has since been used in more than 20 countries across the globe.


By working in partnership with Tetcon, Benmor Medical are able to offer this evacuation equipment to organisations in the UK.

Contact us on 0333 800 9000 to arrange a free demonstration of our evacuation equipment!



For additional resources and to see the evacuation equipment in use please visit:

<https://s-capeplus.com/en/tetcon/>



# S-CAPEPOD

## bariatric evacuation device



The S-CAPEPOD provides a quick and safe evacuation of bed-bound and immobile patients in an emergency – where speed and ease of deployment is critical. The S-CAPEPOD is a simple-to-deploy evacuation device, stored below the bed mattress, which when required can be simply unfolded and fitted around the person laid on the bed mattress.

Features	Benefits
<b>500kg (78st) max load</b>	Able to cater for a range of patient weights
<b>Fast and safe evacuation</b>	You are in control of the evacuation of patients in an emergency
<b>Requires little training</b>	Easy to use
<b>Complies with current standards</b>	Ensure your evacuation plans are in place
<b>Virtually no maintenance</b>	Strong hard wearing and easy to clean

Using the folding flap sections at each side and at the foot-end of the bed to create a sturdy bag in which to envelop the patient and mattress.

Hand-hold straps at the head and foot-end allow the mattress and occupant to be easily maneuvered along the floor and down stairs.

The S-CAPEPOD is suitable for loads up to 500kg (78st).

### Purchase Code

Bariatric S-CAPEPOD 300kg SWL - BMPBSCAPE3  
 Bariatric S-CAPEPOD 500kg SWL - BMPBSCAPE5

### Rental Code

T100SCAPE3  
 T100SCAPE5

Also available for rental with one of our bariatric beds

## technical specification -

<b>Suitable for mattress sizes</b>	Width 100 cm - 125 cm   Max length: 216 cm
<b>Maximum carrying weight</b>	300KG / 500KG
<b>Material</b>	Anti bacterial technology - sanitized technology Fire retardancy - BS 5867/2 Dyneema® fabric is used on the 500kg model
<b>Laundry instructions</b>	60 degrees, dry on low temperature, do not iron, no chemical cleaning, do not bleach, do not use washing powder with fabric softener



By following the correct procedure, the patient is ready for evacuation within 60 seconds

## step by step evacuation -



# S-CAPEPLUS

## emergency evacuation mattress



The S-CAPEPLUS evacuation mattress is designed for the safe evacuation of people that suffer from reduced mobility.

It is multipurpose emergency evacuation support for the physically impaired, such as seniors, wheelchair users, disabled people with injuries, patients with chronological disorders or people with visual impairments.

The S-CAPEPLUS can be used for a horizontal and vertical evacuation. The device is quick to use, user-friendly and comfortable for the patient.

The S-CAPEPLUS is suitable for loads up to 200kg (31.5st).

Features	Benefits
<b>200kg (31.5st) max load</b>	Able to cater for a range of patient weights
<b>Wall Storage</b>	Fast and safe evacuation in a convenient location
<b>Requires very little training</b>	Quick and easy to use
<b>Comfortable and thick foam</b>	Ensuring a safe evacuation without losing comfort
<b>Virtually no maintenance &amp; anti bacterial</b>	Strong hard wearing and easy to clean

Leave no one behind with the S-CAPEPLUS evacuation mattress

**Product code**

BMPBSCAPEPLUS

S-CAPEPLUS 200kg (31.5st) SWL

# technical specification



Ready to use in 30 seconds

<b>Size of mattress</b>	Width 60 cm   Length: 210 cm
<b>Packed size in wall cover</b>	100 cm x 60 cm x 17 cm (L x W x D)
<b>Maximum carrying weight</b>	200KG
<b>Material</b>	Anti bacterial technology - sanitized technology Fire retardancy - BS 5867/2
<b>Laundry instructions</b>	60 degrees, dry on low temperature, do not iron, no chemical cleaning, do not bleach, do not use washing powder with fabric softener



Protected against harmful bacteria



Ergonomic design to avoid discomfort

## step by step evacuation -



- Innovative products & services
- Expert advice
- Highest levels of product conformity
- Design and manufactured products in the UK
- CPD certified training courses
- UK nationwide rental coverage 24/7
- Training and installation
- Patient assessments



#### HEAD OFFICE

The Aurum Centre  
Ham Barn Business Park  
Farnham Road  
Liss  
Hampshire  
GU33 6LB

#### SOUTH WEST BRANCH

The Iridium Centre  
Unit B  
7A4 Victoria Road  
Avonmouth  
BS11 9DB

#### MIDLANDS BRANCH

The Argentum Centre  
Welland Business Park  
Valley Way  
Market Harborough  
Leicestershire  
LE16 7PS

#### NORTHERN BRANCH

The Aereus Centre  
Rainton Bridge Industrial Estate  
Mercantile Road  
Houghton-le-Spring  
Tyne and Wear  
DH4 5PH

supporting bariatric care

Tel: 0333 800 9000  
[www.benmormedical.co.uk](http://www.benmormedical.co.uk)  
[info@benmormedical.co.uk](mailto:info@benmormedical.co.uk)



BS EN 13485:2016 - MD 607484  
CE607483



ISO 9001:2015 - FM 607485