

# Goliath Wheelchair

heavy duty folding bariatric wheelchair



High-strength quality powder-coated steel construction



Detachable leg rests

The Goliath heavy-duty wheelchair is built for strength and durability constructed from high-quality steel tube. The wheelchair is fully folding with a safe working load of up to 300kg (47st).

Unique to the market, the Goliath has the ability of adjusting its seat height, depth and backrest height in order to suit a variety of users. Elevating, height adjustable and detachable arm rests help care givers with patient transfers.

Our wheelchairs come with a rigid, upholstered seat cushion which is fully detachable.

Features	Benefits
<b>300kg (47.3st) max load</b>	Able to cater for a range of patient weights
<b>22" (550mm) to 31.5" (850mm) seat width</b>	Able to cater for a range of patient widths
<b>Folding</b>	Easy storage and transportation
<b>Self propelling wheel</b>	Greater patient independence
<b>Detachable, elevating armrests</b>	Easy side transfers
<b>Detachable leg rests</b>	Allows for patient of varying size
<b>Removable rigid seat cushion</b>	Easy transport, cleaning and storage

Seat Width (mm)	Seat Height (mm)	Seat Depth (mm)	Backrest Height (mm)	Total Width (mm)	Folded Width (mm)	Armrest Height (mm)	Max Load (kg)
550 22"	450-540 18-21.5"	450-540 18-21.5"	420 or 460 17 or 18.5"	750 30"	380 15"	150/250 6/10"	300 47.24st
600 24"	450-540 18-21.5"	450-540 18-21.5"	420 or 460 17 or 18.5"	800 31.5"	380 15"	150/250 6/10"	300 47.24st
650 26"	450-540 18-21.5"	450-540 18-21.5"	420 or 460 17 or 18.5"	850 33.5"	380 15"	150/250 6/10"	300 47.24st
700 28"	450-540 18-21.5"	450-540 18-21.5"	420 or 460 17 or 18.5"	900 35.5"	380 15"	150/250 6/10"	300 47.24st
750 30"	450-540 18-21.5"	450-540 18-21.5"	420 or 460 17 or 18.5"	950 37.5"	380 15"	150/250 6/10"	300 47.24st
800 31.5"	450-540 18-21.5"	450-540 18-21.5"	420 or 460 17 or 18.5"	1000 39.5"	380 15"	150/250 6/10"	300 47.24st
850 31.5"	450-540 18-21.5"	450-540 18-21.5"	420 or 460 17 or 18.5"	1050 41.5"	380 15"	150/250 6/10"	300 47.24st

## Product codes

BMWC008 Goliath folding wheelchair\*

\*A prescription form will be provided for wheelchairs to determine any additional requirements and total price