

Golvo®

patient lift



Key features / spec

200kg (31st) lifting capacity

Quick release hook system - other lifting equipment can be changed quickly and easily

Retractable armrests help to make the perfect mobile lift for gait training during rehabilitation

Vertical lifting motion - adapts automatically to the individual patient and lifting situation

Equilateral and parallel base-width adjustment

Low base - floor clearance of as low as 60mm (2.36")

The Golvo patient lift combines the functionality of an overhead lift with the flexibility of a mobile patient lift.

The parallel base-width adjustment facilitates lifting in many different situations; for example, where under-bed clearance is limited, and lifting to and from chairs.

The Golvo lift offers a very low base and is specially designed for use with beds that have a very low floor clearance.

Parallel and symmetrical base-width adjustment facilitates lifting in many different settings, even where it is otherwise difficult to get close to the patient, while at the same time ensuring optimal stability during the lifting procedure.