

hoist scales

hoist weighing attachment

The hoist weighing attachment has a lightweight design for use on hoists with coat hanger shaped spreader bars, clips on and easy transfer. It is ideal for use within a hospital setting or care home.

A well-padded and comfortable carry case makes the hoist scale genuinely portable from home to home. The clips on the scale make it very easy to move from hoist to hoist to weigh the patient. The hoist scales are now MDD approved and are the only hoist weighing attachments that are legal for use in the UK.



- Lightweight & Portable for use in Home, Community or Hospital
- Suits Most Hoist 'Coathanger' Sling Support
- Safety Clips to Secure

Capacity	200kg
Graduations	100g
Dimensions	630mm x 275mm x 160mm
Weight of scale	6kg
Power Supply	Mains / Rechargeable



Product code: BMW008