

hoist scales

hoist weighing attachment for use on larger hoists

The hoist weighing attachment with a larger design for use on larger hoists with coat hanger shaped spreader bars. The scale clips on and is easy to transfer.

This hoist attachment is designed for use in hospitals and care homes. The clips on the scale make it very easy to move from hoist to hoist to weigh patients.

The hoist scale is now MDD approved which means that they are the only hoist weighing attachment that are legal for use in the UK.



- Lightweight & Portable for use in Home, Community or Hospital
- Suits Most Hoist 'Coathanger' Sling Support
- Safety Clips to Secure
- Fits to Larger Hoists

Capacity	200kg
Graduations	100g
Dimensions	730mm x 275mm x 160mm
Weight of scale	6.7kg
Power Supply	Mains / Rechargeable



Product code: BMW009