

# ReTurn 7600

bariatric sit-to-stand and transfer aid



## Key features / spec

205kg (32st) maximum user weight

Baseplate width 672mm (27"), length 790mm (31")

Complete with wide rising ladder including a central 'power grip'

Transfer platform provides smooth and simple transfers

The height, width and angle of the lower-leg supports can easily be adjusted

Rolls, turns and pivots smoothly and easily on all surfaces

The ReTurn 7600 is a simple to use and incredibly effective assistive device providing convenient and natural support for sit-to-stand and transfer, but it is also a rehabilitation device.

The ReTurn enables safe and active transfers to and from beds, wheelchairs or toilets. Complete with height and width adjustable shin supports that can be used for leg dressings.

For users who still retain some degree of strength, the ReTurn can provide the natural movement during sit-to-stand whilst minimising the risk to the patient or carer.

Optional slings available in various sizes.

### Product codes

BMRT	ReTurn 7600
BMRTBMD	Disposable belt <b>medium</b> (pack of 5)
BMRTBLD	Disposable belt <b>large</b> (pack of 5)
BMRTBXL	Disposable belt <b>extra large</b> (pack of 5)