

rota stand

transfer aid

The Rota Stand has been developed to allow an assisted transfer from one seated position to another. The frame, shin pad and turn table combine to give optimum safety to both carer and patient during the transfer.

The central wheels are a unique innovation, because they allow the Rota Stand to be rolled up to the patient or client even after their feet have been positioned on the product.

BMHE050 rota stand 160kg (25st) weight capacity



rota stand solo

single carer transfer aid

The rotastand solo revolutionises sit to stand and pivot transfers. The T Bar mechanism in the base allows a sit to stand and pivot transfer to be carried out by a single nurse or carer thus reducing the time it takes to carry out the task as well as halving the cost of having two carers/nurses present during the transfer.

The nurse/carers can now assist the client to stand before rotating the client to a new seated position.

BMHE051 rota stand solo 200kg (31st) weight capacity



rota stand compact

foldable transfer aid

The rota stand compact combines all the benefits of a standard rota stand but with the added value of being able to fold yet still lightweight.

The nature of the stand means that it is ideal for community assessments by therapists who often have to lift equipment in and out of their vehicles. Equally, a more independent lifestyle is encouraged for clients who wish to travel to friends and relatives.

The rota stand compact fits easily into the majority of car boots.

BMHE051C rota stand compact 200kg (31st) weight capacity

