

# bariatric commode

standard bariatric commode



## Key features / spec

337kg (53 st) swl

710mm (28") seat width

Includes removable waste container  
with lid

Swing away arms enables  
side transfers

Large space between seat & rear  
handle for wider patients

Plastic bed pan insert (optional)

Our standard bariatric commode can also be used as a shower chair and will fit over a standard sized toilet, making the commode extremely versatile.

Made with an aluminium frame, it has four strong braked castors to ensure stability whilst in use and the swing away arms enable easy side transfers.

**Product code**  
**BMC002**