

ultra twin gantry hoist



Optional UltraTwist

For Smooth Positioning
The UltraTwist accessory offers a unique ability to adjust the patient's position during and after the transfer

Key features / spec

500kg (78st) swl

Range of different rail lengths available

Full 360° rotation of the patient

Freestanding system, not affixed to walls or ceiling

Can be used to lift, transfer and position patients

Can be operated with single hand control using the ultra control unit

The ultra twin freespan gantry hoist can lift, transfer and position patients in a variety of positions. Since it is not affixed to walls or a ceiling, there are no special strength requirements for those structures.

The ultra-twist accessory offers a unique ability to adjust the patient's position during and after the transfer.

Compatible with the aurum bariatric bed.

Product code
BMHUT