

viking xl mobile hoist



optional drop down arms for patient support

Key features / spec

300kg (47st) swl

Weighs only 43kg (7st)

Compatible with the aurum bariatric bed and riser recliner

Wireless hand control included

Equipment height from 49cm - 182cm

Wireless handset included as standard

The viking xl has been developed for bariatric patients but can also be used for smaller patients.

Able to lift patients in a variety of care situations, the viking xl is made of aluminium to be as easy to manoeuvre as possible.

The standard sling bar is the Universal TwinBar 670 which provides more space in the sling for larger patients. The viking can also be equipped with a Liko quick-release hook system to quickly and easily change the sling bar or install a scale.

Product code
BMVXL