

# portable wheelchair weigh beams

300kg weight capacity

The weigh beams are a lightweight construction with a gentle incline on the ramps to enable ease of use. The ends fold in so that the user can easily pack them into a smart, durable carry case. The overall dimensions of the scale accommodate most designs of wheelchair.

The low cost combined with the product features and use in the community makes this the UK's number one selling wheelchair scale.



- Two Portable Weigh Beams
- Lightweight and Portable for use in the Community
- Tare to Deduct Weight of the Wheelchair

Capacity	300kg
Graduations	100g
Dimensions	1200mm x 150mm x 95mm
Weight of scale	8.5kg
Power Supply	Mains / Rechargeable



Product code: BMW003