

X399 bariatric commode

with footplates

This bariatric commode has been designed for comfort and strength and can be wheeled over a standard toilet seat. It is robust and wide enough for bariatric patients, but is also dimensionally sized to fit through a standard door.

The chair has a soft polyurethane padded seat for patient comfort and four 5" braked castors. The commode includes two footplate's that are used to support the patient.



- 254kg (40st) max user weight
- 640mm (25") seat width
- Footplate's clip-on and swing away
- Epoxy powder coated steel
- Swing away armrests
- Easy to clean Polyurethane seat and back
- 585mm (23") seat height - designed to wheel over a standard toilet
- Removable foot complete with footplate

features	benefits
245kg (40st) max weight	Designed to support heavy loads
640mm (25") seat width 460mm (18") seat depth	Wide enough for bariatric patients
Adjustable 'pram handles'	Ensure carer best positioning when transporting patient
Swing away arms	To allow for smooth side transfers
Back framework reinforced	Designed to fit flush with toilet cistern
Seat and Backrest made of Polyurethane	Easy to clean and infection control compliant

Product code
 BMC399