

# BariVest

for assistance with rising,  
standing training and gait training



The BariVest is designed for rising, standing training and gait training with users weighing 180 kg (28st) or more. It can be used in combination with both mobile and stationary lifts, with a four-point sling bar.

BariVest consists of a back section, a front section and a narrower section in the middle to enable a comfortable fit in the groin. It has two sling bar loops, both front and back, as well as belts at three levels.

BariVest is available in two sizes, both for users weighing max 570kg (89st)\*.

\*Maximum user weight of 570kg only applicable to XXXXL BariVest

| features  | benefits   |
|---|--|
| XXXL 570kg (89st) SWL<br>XXXXL 570kg (89st) SWL | Able to be used with heavier patients  |
| Simple design                                   | Easy for the caregiver to apply  |
| Safe for patients and caregivers                | Designed for greater safety and security during raising, standing training and gait training |
| Loop fixing                                     | Can be used with any hook and loop spreader bar  |
| Support in the groin and around the upper body  | Ensuring comfort for the patient   |

| product code | description                          | colour | recommended user weight |
|--------------|--------------------------------------|--------|-------------------------|
| BMBVXXXL     | BariVest, polyester 570kg (89st) swl | grey   | 180-300kg (28-47st)     |
| BMBVXXXXL    | BariVest, polyester 570kg (89st) swl | grey   | 250-570kg (39-89st)     |