

basic sling

bariatric sling for users with good head stability




The basic sling has been developed to provide an effective patient handling solution to the increasing challenges when moving and handling a larger patient.

It has been designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and when the patient has poor torso stability but good head stability.

The sling is designed to be used with any hoist which uses a standard loop attachment spreader bar.

The basic sling is available in two different sizes, xl and xxl and can be used for patients up to 300 kg (47st).

features	benefits
300kg (47st) SWL	Able to be used with heavier patients
Shoulder-high back	The user is supported from the tailbone up to the shoulders and under the legs
Divided leg supports	Easier for the caregiver to apply and remove, the user does not sit on the lifting sling
Loop fixing	Can be used with any hook and loop spreader bar
Upright seated posture	



product code	description	colour	rec user weight	sling seat measurement
BMSXL	basic sling size XL 300kg (47st) swl	grey	100-150 kg (15.5-23.5st)	65cm (26")
BMSXXL	basic sling size XXL 300kg (47st) swl	grey	140-200 kg (22-31.5st)	75cm (30")