

flat sheets



Key features / spec

Handholds down each side of the sheet - ensures carers have a secure grip on the sheet while sliding the patient

Extension strap [optional] eliminates the risk of overstretching and back injury for carers

Different sizes available catering for wider patients.

Developed as an alternative to using tubular slide sheets, the flat sheet incorporates a number of handholds which allow a positive grip when carrying out a patient handling task. The flat sheet can be used singly or in pairs and enjoys the same low friction qualities as the ultra slide sheets.

Product codes

BMHE009	Standard (200cm x 71cm)
BMHE010	Wide (200cm x 90cm)
BMHE011	Handy pack (2 standard flat sheets & 1 compact roller sheet)
BMHE012	Extension strap