

lateral transfer slide

The lateral transfer slide has a smooth surface and a non-slip base allowing a patient to be moved from one surface to another without the need to be lifted. Used in conjunction with the transfer sheet, correct posture during the transfer is maintained. The lateral transfer slide has tapered edges to allow for easy insertion and removal.

Size 185cm x 65cm



transfer slide boards

transfer slide - the low friction surface of the transfer slide allows for a smooth effortless transfer. The advanced construction gives a rigid surface that will withstand weights up to 127kg (20st) without bending and yet is extremely lightweight and easy to pick up and carry. The plastic construction of the transfer slide gives a wipe-clean surface that will not splinter or deteriorate in use.

ultra transfer slide - the ultra transfer slide has the same benefits as the transfer slide, with the addition of a pair of handling loops which give a secure, comfortable grip and allows a patient to pull themselves across the surface of the slide without trapping their fingers underneath.

duo transfer slide - available with or without handles, the unique design of the duo transfer slide allows for an angled transfer as well as a straight transfer from one seated position to another. The handgrips at either end of the deluxe model give the user a means of pulling themselves across the smooth surface while at the same time ensuring the Duo transfer slide is located during the transfer.

- transfer slide 66cm x 20cm and 20st/127kg swl
- ultra transfer slide 66cm x 20cm and 20st/127kg swl
- duo transfer slide



Product codes

BMHE024 lateral transfer slide
 BMHE026 duo transfer slide

BMHE026A
 BMHE027
 BMHE028

duo transfer slide with handles
 transfer slide board
 transfer slide board with handles