

ultra lift pants

bariatric patient movement aid



Key features / spec

Maximum load of 500kg (78st)

Can be used with the Ultra Twin Hoist system (up to 500kg) or the Viking XL hoist (up to 300kg)

All sizes made of net polyester with a pressure-relieving padded crotch

Equipped with two hip straps for an even weight distribution

Available up to size XXXL

Assists with the rehabilitation of a patient

Ultra LiftPants enable caregivers the ability to provide gait training to their patients while providing a sense of security for patients.

Ultra LiftPants provide larger patients with the chance to practice standing and gait training without risk of falling while reducing the risk of injuries to care-givers.

Ultra LiftPants can be used with an overhead or mobile hoist with sufficient lifting height. The advantage of an overhead lift is that the patient can move freely within the lifting area which, in the case of a traverse-type overhead lift, can mean practically the entire room.

Product codes

BMWPXL	XL (500kg)
BMWPXXL	XXL (500kg)
BMWPXXXL	XXXL (500KG)