


repositioning sheet

bariatric repositioning sheet



Key features

500kg (78st) safeworking load

Assists with turning a bariatric patient for personal care

Greatly reduces number of carers needed for repositioning a bariatric patient in bed

Eliminates handling risks to carers

Suitable for use with Aurum Bariatric Bed

Available to hire with mobile or gantry hoist

When a patient is to be moved higher up in the bed or turned, the repositioning sheet is an excellent aid.

The repositioning sheet is designed to be placed under the patient and can be used both in combination with and instead of an ordinary sheet during the entire care period.

A variety of sizes and safe working loads are available.



Product codes

BM3687702

Ultra max 500kg