

transfer sheets



Key features / spec

Long handholds eliminates the need to bend and stretch

Full length sheet provides support for the whole of the patients body

standard transfer sheet

full body-length sheet with six long handholds which can be used with the maxi mover or a bridging board to aid sliding transfers between two surfaces. The long handholds eliminate the need to bend and stretch to effect the transfer. The transfer sheet can also be used in an emergency situation when a patient has fallen on the floor where a hoist is unavailable or inappropriate.

deluxe transfer sheet

The deluxe transfer sheet is used in the same way as the transfer sheet but incorporates a series of extra handholds which allows for more precise positioning of a patient during a transfer. The extra handholds also increase control during the handling task by enabling staff to maintain correct posture during the transfer.

Product codes

BMHE025	Transfer sheet (190cm x 100cm)
BMHE025D	Deluxe transfer sheet (with extra handholds)