

ultra slide sheets



Key features / spec

Large range of sizes available all sizes and shapes of patient

Strong and durable - can be washed and re-used time after time with the same high performance

Can be used for a variety of applications to assist with patient care

Our ultra slide sheets allow three basic patient handling tasks to be carried out without the need to lift and without the risk of sheering to the patient's skin:

- Turning a patient: for x-ray, surgical procedures, recovery, pressure care, changing dressings or cleaning.
- Sitting a patient up in bed: for feeding or relief of chest pain.
- Patient transfers: moving a patient on or off a bed, trolley or examination couch without painful sheering to the hips and shoulders.

Product codes

BMHE001	Standard (122cm x 71cm)
BMHE002	Long (145cm x 71cm)
BMHE003	Wide (122cm x 100cm)
BMHE004	Compact (72cm x 70cm)
BMHE005	Maxi (195cm x 70cm)
BMHE006	Maxi plus (195cm x 140cm)
BMHE007	Mini (60cm x 41cm)
BMHE008	Extra long (200cm x 71cm)