

# quilted slide

## quilted uni slide

Available in two sizes, the uni slide prevents patients from sliding down in a chair or bed. The uni slide allows a one directional slide (back into a chair or up in a bed) then locks, preventing the patient sliding forward. This reduces the amount of patient handling required to reposition the patient and thus reduces the risk of back injury to the nurse or carer associated with these tasks.

The quilted uni slide has a quilted outer surface that gives added comfort for the patient.

**Product codes**

- BMHE190 Quilted Uni-Slide- Chair (45cm x 40cm)
- BMHE191 Quilted Uni-Slide- Bed (80cm x 70cm)



## quilted uni flat slides

Developed as an alternative to the tubular quilted uni slides, the quilted flat uni slide is available in two sizes and incorporates a series of handouts which allow a positive grip when carrying out a patient handling task.

**Product codes**

- BMHE192 Quilted Flat Uni-Slide- Chair (90cm x 40cm)
- BMHE193 Quilted Flat Uni-Slide- Bed (160cm x 70cm)



### Key features

One directional slide	Prevents patient from slipping down in a chair or bed
Different sizes available	Sizes to fit a chair or bed
Quilted slides available	Added comfort for patient