

walking frame

accessories



Key features

Strong construction which can take the weight of patients up to 381kg (60st)

Extremely lightweight and easy to manoeuvre

Available in various sizes to cater for all different types of patient

Various accessories can be fitted to the walking frames to make them better suited to different patients

long leg kit – adds 75mm (3”) to height of frame

wheel kit – 127mm (5”) castors

glide brakes – enables the frame to slide easily across the floor while still being secure when weight is applied to the frame.

leg glides – enables the frame to slide easily across various floor surfaces

All accessories fit to all size frames however fitting any of the accessories reduces the weight capacity of the walking frame. Please contact us for more information.

Product codes

BMWF004	Leg extensions
BMWF005	Wheel kit
BMWF006	Glide brakes
BMWF007	Rear leg glides
BMWF008	Foam Hand Grips (pair)