

# bariatric elbow crutches

suitable for patients up to 318kg (50st)



## Key features

Heavy duty steel

Chromium plated finish

Height adjustable

Non-slip rubber tip

Adjustable from grip to ground and  
armband cuff to handle

Internal cuff size is 14cm

Our heavy duty elbow crutches are suitable for users up to 318kg (50st) and are designed to improve support and comfort for the patient.

The crutches are constructed from heavy duty steel and have a non slip rubber tip.

The crutches are adjustable from grip to ground and armband to cuff.

### Product codes

BMS001AHD

adjustable elbow crutches