

walking stick

single point, adjustable & heavy duty



Key features

Strong steel construction able to take the weight of larger patients

Extremely lightweight and easy to manoeuvre

Height adjustable sticks for general use with different sizes of patient

single point adjustable stick

Height adjustable single point walking stick with a 318kg (50st) weight capacity.

Height adjustable from 29" - 38" to cater for a variety of patients.

Non slip, non-marking rubber tip. Swan neck handle. Suitable for right and left hand use.

Product code

BMS005 Adjustable single point walking stick