

# Phoenix Wheelchair

folding wheelchair



High-strength aluminium frame and double crossbar



Detachable leg rests

The Phoenix wheelchair is built for durability and versatility. Constructed from high-strength aluminum the wheelchair is fully folding with a safe working load of up to 185kg (29st).

Unique to the market, the Phoenix has the ability of adjusting its seat height, depth and backrest height in order to suit a variety of users. Elevating, height adjustable and detachable arm rests help care givers with patient transfers.

Our wheelchairs come with a rigid, upholstered seat cushion with webbing in seat and backrest allowing complete adjustment.

Features	Benefits
185kg (29st) max load	Able to cater for a range of patient weights
23" (590mm) seat width	Able to cater for a range of patient widths
Folding	Easy storage and transportation
Self propelling wheels	Greater patient independence
Detachable, elevating armrests	Easy side transfers
Detachable leg rests	Allows for patient of varying size
Adjustable seat and back rest	Can be used for a variety of patients

Seat Width (mm)	Seat Height (mm)	Seat Depth (mm)	Backrest Height (mm)
590	425-525	420-500	360-460
23.5"	17-21"	17-20"	14.5-18.5"

  

Total Width (mm)	Folded Width (mm)	Armrest Height (mm)	Max Weight (kg)
850	350	140-240	185
33.5"	14"	6/9.5"	29.13st

**Product codes**

T300PW/23

23.5" (590mm) Phoenix folding wheelchair