

bariatric rollator

adjustable bariatric rollator



Optional forearm gutter supports and hand grips



Optional IV pole and oxygen holder

Key features / spec

453kg (71st) swl

Strong steel construction with a durable and reliable frame

Folds for storage and easy transportation

Width adjustable from 660mm to 812mm (26" - 32")

Height adjustable from 762mm - 889mm (30" - 35")

Pressure activated braking system

Our bariatric rollator offers stability, safety and freedom of movement for the larger patient and can meet the needs of people who have difficulty with their balance.

The rollator is constructed from steel with strong welded joints. It is lightweight and easy to store. The castors are skid resistant and easy rolling with a braking system on the rear wheels which is activated by applying weight to the frame.

A range of accessories are also available including detachable seats, baskets and gutter supports.

Product codes

BMR001	bariatric rollator
BMR002	forearm gutter supports
BMR004	IV pole
BMR005	oxygen holder