

Goliath Wheelchair

heavy duty folding bariatric wheelchair



High-strength quality powder-coated steel construction



Detachable leg rests

The Goliath heavy-duty wheelchair is built for strength and durability constructed from high-quality steel tube. The wheelchair is fully folding with a safe working load of up to 300kg (47st).

Unique to the market, the Goliath has the ability of adjusting its seat height, depth and backrest height in order to suit a variety of users. Elevating, height adjustable and detachable arm rests help care givers with patient transfers.

Our wheelchairs come with a rigid, upholstered seat cushion which is fully detachable.

Features	Benefits
300kg (47.3st) max load	Able to cater for a range of patient weights
27.5" (700mm) or 29.5" (750mm) seat width	Able to cater for a range of patient widths
Folding	Easy storage and transportation
Self propelling wheels	Greater patient independence
Detachable, elevating armrests	Easy side transfers
Detachable leg rests	Allows for patient of varying size
Removable rigid seat cushion	Easy transport, cleaning and storage

Seat Width (mm)	Seat Height (mm)	Seat Depth (mm)	Backrest Height (mm)
700	450-540	450-540	420 or 460
28"	18-21.5"	18-21.5"	17 or 18.5
750	450-540	450-540	420 or 460
30"	18-21.5"	18-21.5"	17 or 18.5"
Total Width (mm)	Folded Width (mm)	Armrest Height (mm)	Max Weight (kg)
900	380	150-250	300
35.5"	15"	6/10"	47.24st
950	380	150-250	300
37.5"	15"	6/10"	47.24"

Product codes

- T500GW28 28" (700mm) Goliath heavy duty folding wheelchair
- T500GW30 30" (750mm) Goliath heavy duty folding wheelchair
- T500GW31 31" (787mm) Goliath heavy duty folding wheelchair

supporting bariatric care