

Freeway Commode / Shower Chair



Padded seat for optimum comfort



Adjustable heel supports which are removable and swing away

The Freeway can be used as both a commode and a shower chair. Constructed from increased diameter tubing and heavy duty castors, this commode can support the weight of bariatric users up to 350kg (55st).

This commode/ shower chair can be made in a variety of seat widths ranging from 590mm (23") to 890mm (35") to ensure best fit for patients.

With the added benefit of padded seats and height adjustable heel supports, the Freeway Commode is comfortable for users.

Features	Benefits
350kg (55st) Safe Working Load	Caters for bariatric patients
Range of seat widths	Ensures best fit for the user
Commode/ shower chair	Suitable for a variety of applications for bariatric patients
FREECOTE antimicrobial and antibacterial coating	Resistant against MRSA, e-coli, salmonella and listeria.

Product codes

BMCT100590 - Freeway Commode 590mm (23") seat width
BMCT100640 - Freeway Commode 640mm (25") seat width
BMCT100690 - Freeway Commode 690mm (27") seat width
BMCT100740 - Freeway Commode 740mm (29") seat width
BMCT100840 - Freeway Commode 840mm (33") seat width
BMCT100890 - Freeway Commode 890mm (35") seat width
BMCT100CP - Commode Pan for Freeway