

Tilt In Space Riser Recliner

available for rental



25 stone
SWL



Key features / spec

159kg (25st) swl

Seat width 19"

Waterfall back

Waterproof and infection resistant vinyl fabric for easy cleaning

Mobile hoist compatible

Dual motor allows for individual electric control of all functions

To promote better posture for the user, the tilt in space function ensures the angle of the seat to back remains the same whilst the chair is raised or reclined

The chair is mounted on 4 heavy duty castors and includes handles on the rear of the chair for easy transportation.

Space underneath the chair can accommodate mobile hoists for easier patient transfers.

Rental code
T100TIS/25