

Flojac

air assisted lifting device

500kg
SWL



Key features / spec

500kg SWL - Can accommodate a large range of patients

Head end support wedge - Provides added comfort for the patient

Multi layered - For added pressure relief

Wipeable - Meets infection control standards

Simple pump design interface - Easy to inflate and deflate

The Flojac is an air assisted lifting device, designed to help clinicians safely lift patients weighing up to 500kg from the floor.

The powered device can also reduce the manual handling required for lateral transfers from one surface to another, reducing the risk of musculoskeletal injury to staff.

The pump is operated by a simple on/off switch, making it quick and easy to inflate and deflate. It is equipped with a head support wedge for added patient comfort and 10 TPU, wipeable handles to encourage better posture for caregivers.

Rental code
T100FLO