

# folding wheelchair scales

with 300kg weight capacity

The folding wheelchair scale has a large weighing platform to accommodate almost every wheelchair and a gentle incline on the ramps to enable ease of use.

The ramps fold in and the platform folds up and onto a trolley so that the wheelchair scales can be easily and safely transported, as well as providing a compact unit for easy storage. This scale has a tare and pre-set tare function so that the weight of the wheelchair can be deducted before weighing the patient.



- Low Profile Platform
- Tare to Deduct Weight of Wheelchair
- Folded for Easy Transport and Storage
- Ideal for Carehomes



Capacity	300kg
Graduations	100g
Dimensions	700mm x 800mm x 80mm
Weight of scale	38kg
Power Supply	Mains / Rechargeable



Product code: BMW001