


# pro-lock slings

## amputee disposable sling

The pro-lock amputee disposable sling has been developed to provide a low cost effective patient handling solution to the increasing challenges of infection control. Available in a range of sizes, the pro-lock amputee sling is designed to be issued for single patient use only and to be disposed of when soiled or no longer required.

The pro-lock amputee sling incorporates a head support for increased comfort and a patented positive locking clip fixing, suitable for use with any hoist which uses a clip fixing mechanism.



### features

Patient specific	Infection control maximised
Unique safety harness with positive locking clip fixing	Greater security for the patient
Low cost	Laundry costs eliminated
230kg/36st safe working load	Able to be used with heavier patients

product code	size	dimensions
BMS304XS	extra small	230kg (36st)
BMS304	small	230kg (36st)
BMS305	medium	230kg (36st)
BMS306	large	230kg (36st)
BMS306XL	extra large	230kg (36st)