


# pro-lock

## patient specific slings


The pro-lock disposable sling has been developed to provide a low cost effective patient handling solution to the increasing challenges of infection control. Available in a range of sizes, the pro-lock is designed to be issued for single patient use only and to be disposed of when soiled or no longer required.

The pro-lock sling incorporates a head support for increased comfort and a patented positive locking clip fixing, suitable for use with any hoist which uses a clip fixing mechanism.



### features

Patient specific	Infection control maximised
Unique safety harness with positive locking clip fixing	Greater security for the patient
Low cost	Laundry costs eliminated
230kg/36st safe working load	Able to be used with heavier patients



product code	size	strap colour
BMS195	extra small	brown
BMS196	small	red
BMS197	medium	yellow
BMS198	large	green
BMS198XL	extra large	grey

Also available as pro-lock amputee slings