

Civitas² Bariatric Community Bed

Electrically operated backrest, height,
tilt and knee break

380kg
SWL



Key features / spec

380kg safe working load, with a maximum user weight of 350kg (55 stone)

Electrically operated backrest, height, tilt and knee break

Integral extending platform - catering for taller patients

Auto contour and auto regression functions

Low height of 30cm can be achieved

Compatible with a variety of mattresses

With an attractive domestic appearance, the civitas² bariatric community bed can accommodate a variety of patients weighing up to 350kg.

With a low mattress platform height ranging from 30cm to 70cm and an integral extending platform of 15cm, the community bed is extremely versatile.

The fully profiling bed has an easy to use handset featuring a lockout function for patient safety. It features drop down, detachable side rails for manual handling activities.

Product code
BMAPH188